en Sterilisation Guidance for Procedure Pack manufacturers

This product is intended for sale to manufacturers of procedure packs. It is offered NON-STERILE in bulk packaging and must be sterilised prior to the distribution as a finished product. The procedure pack manufacturer is responsible to validate and verify their chosen packaging and sterilisation conditions to ensure the finished product is sterile. The procedure pack manufacturer is responsible to label the finished product with the sterilisation method used. The recommended sterilisation methods are Ethylene Oxide (EO). When Ethylene Oxide is used, the process shall comply with ISO 11135 and ISO 10993-7. Product performance has been evaluated for the following conditions: **Ethylene oxide:** Maximum cycle time: 168 hours, Humidity (RH): 45–90%, Maximum temperature: 55°C, Minimum pressure: 60 mBarA. Mölnlycke provides a non-sterile product intended to be sterilised before use. If sterilised in conditions outside what is stated above, the product characteristics and performance may be adversely affected in such a way that the health and safety of the patient or user might be compromised.

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