BORT RhizoFlex Daumen-Ringorthese

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REF 105 100

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D105100|2023-07|007 ML | Rev. 01

DE BORT RhizoFlex Daumen-Ringorthese

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EN BORT RhizoFlex Thumb ring brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose
This medical device is a thumb brace for immobilisation of the thumb saddle joint (CMC joint).

Indications
Post-operative, post-traumatic, or degenerative irritations in the thumb saddle joint area, e.g. arthrosis of the CMC-I joint (thumb saddle joint), post-operatively following arthroplasty of the CMC-I joint, instability of the CMC-I joint.

Contraindications
Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes
This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
– remove the medical device prior to radiological examinations
– do not close it too tightly during longer periods of rest, e.g. when sleeping
– in case of numbness, loosen the medical device and remove it if necessary
– in case of persistent complaints, consult the physician or specialist retailer
– use the medical device in accordance with therapeutic needs
– only use other products simultaneously after consultation with your physician
– do not make any changes to the product
– do not wear it on open wounds
– do not use in case of intolerance of one of the materials used
– no re-use. This medical aid is intended for treating one patient only
– whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions
1. Slide the brace over your thumb, the guide loop points upwards here. The aluminium stabilisation element integrated in the brace frame incorporates the base of the thumb. The Velcro fastener is open here, however it remains in the guide loop.
2. Then close the Velcro tape.
3. As the aluminium stabilisation element is formable, this brace can repeatedly be refitted during the course of therapy in order to retain an optimum fit.

Removal
To remove the brace, please open the circumferential strap and slide the brace off over the thumb.

Material composition
Plastic frame: 70% thermoplastic polyurethane (TPU), 30% thermoplastic elastomers (TPE)
Sealing tape: 55% polyester, (PES), 45% polyamide (PA)

Cleaning information
 Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. If required, clean the brace with soap and lukewarm water. Air dry the brace.

Guarantee
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

MD Medical device | Single patient – multiple use

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