
Mölnlycke®

Mepilex® Transfer Ag

WITH SAFETAC® TECHNOLOGY

STERILE | **EO**

C **E** 2797



UK
CA 0086

Manufacturer / Fabricant

 **Mölnlycke Health Care AB**
Entreprenörstråket 21
SE-431 53 Mölndal, Sweden

www.molnlycke.com

Mepilex® Transfer Ag

www.molnlycke.com/symbols

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Antimicrobial soft silicone exudate transfer dressing

Product description

Mepilex Transfer Ag consists of a Safetac® wound contact layer, and a compressed polyurethane foam containing silver sulphate and coloured by activated carbon.

Safetac® is a proprietary adhesive technology developed by Mölnlycke that minimises pain to patients and trauma to intact skin or wounds.

Mepilex Transfer Ag is a soft silicone wound contact layer that absorbs and transfers exudate, maintains a moist wound healing environment and has antimicrobial properties.

Mepilex Transfer Ag is soft and conformable making it easy to keep the dressing in contact with the wound surface and the surrounding skin, even in awkward or uneven areas. Mepilex Transfer Ag maintains a moist wound environment in combination with an appropriate secondary dressing.

Mepilex Transfer Ag contains silver sulphate to create an effective bacterial barrier and inactivate a wide range of wound related pathogens (bacteria and fungi), as shown *in vitro*. By reducing the number of micro-organisms, Mepilex Transfer Ag may also reduce odour.

Dressing material content:
Polyurethane, silicone, polyolefin, polyacrylate, silver sulphate which equals to 1.2 mg/cm² silver

Indications for use

Mepilex Transfer Ag is designed for a wide range of low to high exuding leg and foot ulcers, pressure ulcers, malignant wounds, partial thickness burns and donor sites, where there is a risk of infection.

Mepilex Transfer Ag may be used on infected wounds as part of a treatment regimen under supervision of a qualified health care professional.

Mepilex Transfer Ag can be used under compression bandaging.

Precaution(s)

- Do not use Mepilex Transfer Ag on patient and/or user with known hypersensitivity to the ingoing materials/components of the product.
- Clinicians/healthcare professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Frequent or prolonged use of this product may result in transient or permanent discoloration of skin.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

Instructions for use

Mepilex Transfer Ag should be used by or under the supervision of a qualified health care professional.

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Select appropriate dressing size. For best clinical outcomes, ensure that the product covers the wound bed and overlaps the wound edges. If more than one piece of Mepilex Transfer Ag is required overlap the edges.
4. Cut to appropriate size (if needed). Remove the release films and apply the product with the adherent side to the wound. Do not stretch.
5. Apply a secondary dressing. Choose secondary dressing appropriate for the level of exudate. The secondary dressing should exceed the edges of Mepilex Transfer Ag. Fixate to secure close contact with the wound.

Each Mepilex Transfer Ag dressing may be left in place for up to 14 days depending on the patient, condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

Mepilex Transfer Ag is intended for short-term use up to 4 weeks. For continued use, a re-assessment by a physician is recommended. In case of clinical infection Mepilex Transfer Ag does not replace the need for systemic or other adequate infection treatment.

The interaction of Mepilex Transfer Ag together with cleansing agents, oxidising agents such as hypochlorite solutions/hydrogen peroxide and topical treatments has not been demonstrated.

Prior to commencing radiation therapy, remove Mepilex Transfer Ag if product is present in the treatment area. A new dressing can be applied following treatment.

Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).

Special storage conditions and handling conditions

Mepilex Transfer Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any colour variances in Mepilex Transfer Ag do not affect product performance or safety.

Disposal should be handled according to local environmental procedures.

Other information

For EU-customers, further information on the safety and clinical performance of Mepilex Transfer Ag please refer to the 'Summary of Safety and Clinical Performance' for Mepilex Transfer Ag. The document is published in the European Database on Medical Devices, EUDAMED, <https://ec.europa.eu/tools/eudamed>.

If any serious incident has occurred in relation to the use of Mepilex Transfer Ag, you should report to Mölnlycke Health Care and to your local competent authority.

Mepilex and Safetac are registered trademarks of Mölnlycke Health Care AB.

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