













Orbit[®]soft

Infusion Set

en	Instructions for use	10
ar	إرشادات الاستخدام	17
cs	Návod k použití	18
da	Brugsanvisning	22
de	Gebrauchsanweisung	26
el	Οδηγίες χρήσης	30
es	Instrucciones de uso	34
et	Kasutusjuhend	38
fi	Käyttöohje	42
fr	Mode d'emploi	46
he	הוראות שימוש	53
hr	Upute za uporabu	54
hu	Használati útmutató	58
it	Istruzioni per l'uso	62
ja	取扱説明書	66
ko	사용 설명서	70
lt	Naudojimo instrukcija	74
lv	Lietošanas instrukcija	78
nl	Gebruiksaanwijzing	82
no	Bruksanvisning	86
pl	Instrukcja obsługi	90
pt	Instruções de utilização	94
ro	Instrucțiuni de utilizare	98
ru	Инструкция по использованию	102
sk	Návod na použitie	106
sl	Navodila za uporabo	110
sv	Bruksanvisning	114
tr	Kullanım talimatları	118
zh	使用说明	122

STERILE EO

en	Sterilized using ethylene oxide
ar	مُعقمة باستخدام أكسيد الإيثيلين
cs	Sterilizováno ethylenoxidem
da	Steriliseret med etylenoxid
de	Sterilisiert mit Ethylenoxid
el	Αποστειρωμένο με αιθυλενοξειδίο
es	Esterilizado con óxido de etileno
et	Steriliseeritud etüleenoksiidiga
fi	Steriloitu etyleenioksidilla
fr	Stérilisé avec de l'oxyde d'éthylène
he	מעוקר באמצעות אתילן-אוקסיד
hr	Sterilizirano etilen oksidom
hu	Etilén-oxiddal sterilizálva
it	Sterilizzato mediante ossido di etilene
ja	酸化エチレンを用いて滅菌済み
ko	에틸렌옥사이드로 소독
lt	Sterilizuota naudojant etileno oksidą
lv	Sterilizēts ar etilēna oksīdu
nl	Gesteriliseerd met behulp van ethyleenoxide
no	Sterilisert med etylenoksid
pl	Wysterylizowany tlenkiem etylenu
pt	Esterilizado com óxido de etileno
ro	Sterilizat cu etilenoxid
ru	Стерилизовано этиленоксидом
sk	Sterilizované etylénoxidom
sl	Sterilizirano z etilenoksidom
sv	Steriliserad med etylenoxid
tr	Etilen oksit ile sterilize edilmiştir
zh	使用环氧乙烷消毒

en	Single sterile barrier system. Sterilized using ethylene oxide.
ar	نظام الحاجز الواحد المعقم. معقم باستخدام أكسيد الإيثيلين.
cs	Systém jednovrstvé sterilní bariéry. Sterilizováno ethylenoxidem.
da	Enkelt steril barrieresystem. Steriliseret med etylenoxid.
de	Sterilbarrieresystem. Mit Ethylenoxid sterilisiert.
el	Σύστημα μονού στείρου φραγμού. Αποστειρωμένο με αιθυλενοξειδίο.
es	Sistema de barrera estéril simple. Esterilizado con óxido de etileno.
et	Ühekordne steriilse barjääri süsteem. Steriliseeritud etüleenoksiidiga.
fi	Yksittäinen steriili estojärjestelmä. Steriloitu etyleenioksidilla.
fr	Système de barrière stérile simple. Stérilisé avec de l'oxyde d'éthylène.
he	מערכת חסם סטריילי יחיד. מעוקר באמצעות אתילן-אוקסיד.
hr	Sustav s jednostrukom sterilnom barijerom. Sterilizirano etilen oksidom.
hu	Egyszeres steril gátrendeszer. Etilén-oxiddal sterilizálva.
it	Sistema di barriera sterile singolo. Sterilizzato mediante ossido di etilene.
ja	単独使用無菌バリアシステム。酸化エチレン滅菌処理済み。
ko	단일 멸균 차단 시스템. 에틸렌 옥사이드를 사용하여 멸균합니다.
lt	Vienguba sterili barjero sistema. Sterilizuota naudojant etileno oksidą.
lv	Vienas sterilas barjeras sistēma. Sterilizēta ar etilēna oksīdu.
nl	Enkel steriel barrieresysteem. Gesteriliseerd met behulp van ethyleenoxide.
no	System med én steril barriere. Sterilisert med etylenoksid.
pl	System pojedynczej bariery sterylnej. Wysterylizowany tlenkiem etylenu.
pt	Sistema de barreira estéril simples. Esterilizado com óxido de etileno.
ro	Sistem cu barieră sterilă unică. Sterilizat cu etilenoxid.
ru	Единая барьерная система стерилизации. Стерилизовано этиленоксидом.
sk	Jednoduchý sterilný bariérový systém. Sterilizované etylénoxidom.
sl	Enojni sterilni pregradni sistem. Sterilizirano z etilenoksidom.
sv	Enkelt sterilt barriärsystem. Steriliserad med etylenoxid.
tr	Tek steril bariyer sistemi. Etilen oksit ile sterilize edilmiştir.
zh	单次无菌屏障系统。使用环氧乙烷灭菌。

en	Consult instructions for use
ar	اتبع إرشادات الاستخدام
cs	Prostudujte si návod k použití
da	Se venligst brugervejledningen
de	Gebrauchsanweisung beachten
el	Ανατρέξτε στις οδηγίες χρήσης
es	Consulte las instrucciones de uso
et	Järgige kasutusjuhendit
fi	Katso käyttöohjeita
fr	Consultez les instructions d'utilisation
he	יש לעיין בהוראות השימוש
hr	Pogledajte upute za uporabu
hu	Lásd a használati utasítást
it	Consultare le istruzioni per l'uso
ja	取扱説明書を参照のこと
ko	사용 설명서를 참고할 것
lt	Skaitykite naudojimo instrukciją
lv	levērojiet lietošanas instrukciju
nl	Raadpleeg gebruiksaanwijzing
no	Følg bruksanvisningen
pl	Przestrzegać instrukcji obsługi
pt	Consultar as instruções de utilização
ro	A se consulta instrucțiunile de utilizare
ru	Обратитесь к инструкции по применению
sk	Pozrite si návod na používanie
sl	Upoštevajte navodila za uporabo
sv	Följ bruksanvisningen
tr	Kullanım talimatlarını dikkate alın
zh	参见使用说明

en	Do not re-use
ar	لا تقم بإعادة الاستخدام مرة أخرى
cs	Nepoužívejte opakovaně
da	Må ikke genanvendes
de	Nicht wiederverwenden
el	Μην επαναχρησιμοποιείτε
es	No reutilizar
et	Ärge kasutage korduvalt
fi	Ei saa käyttää uudelleen
fr	Ne pas réutiliser
he	לא לשימוש חוזר
hr	Nije za višekratnu uporabu
hu	Ne használja újra
it	Non riutilizzare
ja	再使用しないでください
ko	재사용하지 마십시오
lt	Negalima naudoti pakartotinai
lv	Nelietot atkārtoti
nl	Niet hergebruiken
no	Kun for engangsbruk
pl	Nie używać ponownie
pt	Não reutilizar
ro	A nu se reutiliza
ru	Не использовать повторно
sk	Nepoužívajte opakovane
sl	Ni za ponovno uporabo
sv	Endast för engångsbruk
tr	Tekrar kullanmayın
zh	不得重复使用

en	Non-pyrogenic
ar	لا تؤدي للإصابة بالحمى
cs	Nepyrogeenní
da	Non-pyrogen
de	Pyrogenfrei
el	Μη πυρετογόνο
es	Apirógeno
et	Mittepürogeenne
fi	Pyrogeeniton
fr	Apyrogène
he	לא פירוגני
hr	Nepirogeno
hu	Nem pirogén
it	Apirogeno
ja	発熱性物質を含まない
ko	비발열성
lt	Nepirogeniškas
lv	Apirogēns
nl	Niet-pyrogeen
no	Ikke-pyrogen
pl	Niepirogenny
pt	Não-piragénico
ro	Apirogen
ru	Апирогенно
sk	Nepyrogeenne
sl	Apirogeno
sv	Pyrogenfri
tr	Pirojenik değildir
zh	无热原

en	Do not use if package is damaged
ar	لا تستخدم المنتج في حالة تلف العبوة
cs	Nepoužívejte, pokud je obal poškozený
da	Må ikke anvendes, hvis emballagen er beskadiget
de	Nicht verwenden, wenn Verpackung beschädigt ist
el	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιού
es	No utilizar si el embalaje está dañado
et	Ärge kasutage, kui pakend on kahjustatud
fi	Ei saa käyttää, jos pakkaus on vaurioitunut
fr	Ne pas utiliser si l'emballage est endommagé
he	אין להשתמש אם האריזה ניזוקה
hr	Ne koristiti ako je pakiranje oštećeno
hu	Ne használja, ha a csomagolás megsérült.
it	Non utilizzare se la confezione è danneggiata
ja	包装に損傷がある場合は、使用しないでください
ko	패키지가 손상된 경우 사용하지 마십시오
lt	Nenaudoti, jei pakuotė pažeista
lv	Nelietot, ja iepakojums ir bojāts
nl	Niet gebruiken als verpakking beschadigd is
no	Skal ikke brukes hvis emballasjen er skadet
pl	Nie używać, jeśli opakowanie jest uszkodzone
pt	Não utilizar se a embalagem estiver danificada
ro	A nu se utiliza dacă ambalajul este deteriorat
ru	Не использовать при повреждении упаковки
sk	Nepoužíte, ak je obal poškodený.
sl	Ne uporabljajte, če je ovojnina poškodovana
sv	Använd inte produkten om förpackningen är skadad
tr	Ambalajı hasarlıysa kullanmayın
zh	若包装破损不得使用

MD

en	Medical device
ar	جهاز طبي
cs	Zdravotnický prostředek
da	Medicinsk udstyr
de	Medizinprodukt
el	Ιατροτεχνολογικό προϊόν
es	Producto sanitario
et	Meditsiiniline seade
fi	Lääkinnällinen laite
fr	Dispositif médical
he	מכשור רפואי
hr	Medicinski proizvod
hu	Gyógyászati eszköz
it	Dispositivo medico
ja	医療機器
ko	의료 기기
lt	Medicinos prietaisas
lv	Medicīniskā ierīce
nl	Medisch hulpmiddel
no	Medisinsk utstyr
pl	Urządzenie medyczne
pt	Dispositivo médico
ro	Dispozitiv medical
ru	Медицинское изделие
sk	Zdravotnícka pomôcka
sl	Medicinski pripomoček
sv	Medicinteknisk produkt
tr	Tıbbi cihaz
zh	医疗器械

en	Caution
ar	تحذير
cs	Upozornění
da	Forsigtig
de	Achtung
el	Προσοχή
es	Precaución
et	Ettevaatust!
fi	Huomio
fr	Attention
he	אזהרה
hr	Oprez
hu	Figyelem
it	Attenzione
ja	注意
ko	주의
lt	Atsargiai
lv	Uzmanību!
nl	Voorzichtig
no	Forsiktig
pl	Przestroga
pt	Cuidado
ro	Atenție
ru	Осторожно!
sk	Upozornenie
sl	Svarilo
sv	Försiktighet
tr	Dikkat
zh	注意

en	Contains phthalates
ar	يحتوي على فثالات
cs	Obsahuje ftaláty
da	Indeholder ftalater
de	Enthält Phthalate
el	Περιέχει φθαλικές ενώσεις
es	Contiene ftalatos
et	Sisaldab ftalaate
fi	Sisältää ftalaatteja
fr	Contient des phthalates
he	מכיל פתאלאט
hr	Sadržava ftalate
hu	Ftalátot tartalmaz
it	Contiene ftalati
ja	フタル酸エステル類を含有
ko	프탈레이트 포함
lt	Sudėtyje yra ftalatų
lv	Satur ftalātus
nl	Bevat ftalaten
no	Inneholder ftalater
pl	Zawiera ftalany
pt	Contém ftalatos
ro	Conține ftalați
ru	Содержит фталаты
sk	Obsahuje ftaláty
sl	Vsebuje ftalate
sv	Innehåller ftalater
tr	Ftalat içerir
zh	含邻苯二甲酸盐

REF

en	Catalogue number
ar	الرقم الكتالوجي
cs	Katalogové číslo
da	Art.nr.
de	Artikelnummer
el	Αρ. καταλόγου
es	Número de catálogo
et	Tellimisnumber
fi	Tuotenumero
fr	Référence catalogue
he	מספר קטלוגי
hr	Kataloški broj
hu	Megrendelési szám
it	Riferimento di catalogo
ja	カタログ番号
ko	주문 번호
lt	Užsakymo numeris
lv	Pasūtījuma numurs
nl	Artikelnummer
no	Artikkelnummer
pl	Numer katalogowy
pt	Número de catálogo
ro	Număr de catalog
ru	Номер по каталогу
sk	Katalógové číslo
sl	Kataloška številka
sv	Referensnummer
tr	Katalog numarası
zh	目录编号

**LOT**

en	Batch Code
ar	رمز الدفعة
cs	Číslo šarže
da	Lotnummer
de	Fertigungslosnummer
el	Κωδικός παρτίδας
es	Código de lote
et	Partii number
fi	Eränumero
fr	Numéro de lot
he	קוד אצווה
hr	Broj serije
hu	Gyártási tételszám
it	Codice di lotto
ja	ロット番号
ko	배치 번호
lt	Partijos kodas
lv	Partijas numurs
nl	Batch-code
no	Batchkode
pl	Numer partii
pt	Código do lote
ro	Codul lotului
ru	Код партии
sk	Číslo šarže
sl	Številka serije
sv	Lotnummer
tr	Parti Kodu
zh	批号

en	Date of manufacture
ar	تاريخ التصنيع
cs	Datum výroby
da	Fremstillingsdato
de	Herstellungsdatum
el	Ημερομηνία κατασκευής
es	Fecha de fabricación
et	Tootmiskuupäev
fi	Valmistuspäivämäärä
fr	Date de fabrication
he	תאריך ייצור
hr	Datum proizvodnje
hu	Gyártás ideje
it	Data di fabbricazione
ja	製造日
ko	제조일
lt	Pagaminimo data
lv	Izgatavošanas datums
nl	Productiedatum
no	Produksjonsdato
pl	Data produkcji
pt	Data de fabrico
ro	Data fabricației
ru	Дата изготовления
sk	Deň výroby
sl	Datum izdelave
sv	Tillverkningsdatum
tr	Üretim tarihi
zh	制造日期

en	Use by date
ar	تاريخ الصلاحية
cs	Datum použitelnosti
da	Udløbsdato
de	Verwendbar bis
el	Ημερομηνία λήξης
es	Fecha de caducidad
et	Kõlblik kuni
fi	Viimeinen käyttöpäivämäärä
fr	Date limite d'utilisation
he	תאריך תפוגה
hr	Upotrijebiti do
hu	Felhasználható
it	Data di scadenza
ja	使用期限
ko	사용기한
lt	Tinka naudoti iki
lv	Derīgs līdz
nl	Uiterste gebruiksdatum
no	Utløpsdato
pl	Termin przydatności
pt	Data de validade
ro	Data de expirare
ru	Использовать до
sk	Použiteľné do
sl	Uporabno do
sv	Utgångsdatum
tr	Son kullanım tarihi
zh	使用截止日期

en	Manufacturer
ar	الجهة المصنعة
cs	Výrobce
da	Producent
de	Hersteller
el	Κατασκευαστής
es	Fabricante
et	Tootja
fi	Valmistaja
fr	Fabricant
he	יצרן
hr	Proizvođač
hu	Gyártó
it	Fabbricante
ja	製造元
ko	제조업체
lt	Gamintojas
lv	Ražotājs
nl	Fabrikant
no	Produsent
pl	Producent
pt	Fabricante
ro	Producător
ru	Изготовитель
sk	Výrobca
sl	Izdelovalec
sv	Tillverkare
tr	Üretici Firma
zh	生产商

GTIN**UDI**

en	Unique Device Identifier	en	Keep dry
ar	المعرف المميز للجهاز	ar	احفظ المنتج جافاً
cs	Jedinečný identifikátor prostředku	cs	Uchovávejte v suchu
da	Unik udstyrsidentifikation	da	Opbevares tørt
de	Produktidentifizierungsnummer	de	Trocken aufbewahren
el	Μοναδικό αναγνωριστικό προϊόντος	el	Διατηρείτε στεγνό
es	Identificador único del dispositivo	es	Mantener seco
et	Unikaalne seadme identifikaator	et	Hoida kuivas kohas
fi	Yksilöllinen laitetunniste	fi	Pidettävä kuivana
fr	Identifiant unique des dispositifs médicaux	fr	Conserver au sec
he	מזהה מכשיר ייחודי	he	שמור יבש
hr	Jedinstveni broj trgovačke jedinice	hr	Čuvati na suhom mjestu
hu	Egyedi eszközazonosító	hu	Tartsa szárazon
it	Identificativo univoco del dispositivo	it	Conservare in un luogo asciutto
ja	UDI情報	ja	雨のあたらない場所に保管すること
ko	고유한 기기 식별 번호	ko	빗물을 멀리할 것
lt	Unikalus įrenginio identifikatorius	lt	Laikyti sausoje vietoje
lv	Unikāls ierīces identifikators	lv	Uzglabāt sausā vietā
nl	Unieke apparaatidentificatie	nl	Droog bewaren
no	Unik enhetsidentifikator	no	Oppbevares tørt
pl	Unikalny identyfikator urządzenia	pl	Chronić przed wilgocią
pt	Identificação única do dispositivo	pt	Conservar seco
ro	Identificator unic al dispozitivului	ro	A se păstra uscat
ru	Уникальный идентификатор изделия	ru	Беречь от влаги
sk	Jedinečný identifikátor zariadenia	sk	Uchovávať v suchu
sl	Edinstveni identifikator pripomočka	sl	Hranite na suhem
sv	Unik enhetsidentifikation	sv	Förvaras torrt
tr	Benzersiz cihaz tanımlayıcı	tr	Kuru yerde tutun
zh	唯一设备标识符	zh	保持干燥

en	Keep away from sunlight	en	Distributor
ar	يحفظ بعيداً عن ضوء الشمس	ar	الموزع
cs	Chraňte před slunečním zářením	cs	Distributor
da	Må ikke udsættes for direkte sollys	da	Distributor
de	Vor Sonnenlicht schützen	de	Distributor
el	Φυλάσσετε μακριά από το ηλιακό φως	el	Διανομέας
es	Mantener alejado de la luz solar	es	Distribuidor
et	Hoida eemal päikesevalgusest	et	Edasimüüja
fi	Pidettävä poissa auringonvalosta	fi	Jälleenmyyjä
fr	Conserver à l'abri de la lumière du soleil	fr	Distributeur
he	שמור הרחק מאור השמש	he	מפיץ
hr	Držati zaštićeno od sunčeve svjetlosti	hr	Distributer
hu	Tartsa napfénytől távol	hu	Forgalmazó
it	Conservare al riparo dalla luce solare	it	Distributore
ja	直射日光を避けて保管すること	ja	販売元
ko	직사광선을 피할 것	ko	공급자
lt	Saugoti nuo saulės spindulių	lt	Platintojas
lv	Sargāt no saules gaismas	lv	Izplatītājs
nl	Niet blootstellen aan zonlicht	nl	Distributeur
no	Beskyttes mot sollys	no	Distributør
pl	Chronić przed światłem słonecznym	pl	Dystrybutor
pt	Manter afastado da luz solar	pt	Distribuidor
ro	A se păstra departe de acțiunea soarelui	ro	Distribuitor
ru	Беречь от солнечных лучей	ru	Дистрибутор
sk	Chránite pred slnečným svetlom	sk	Distribútor
sl	Hranite zaščiteno pred sončno svetlobo	sl	Distributer
sv	Skyddas från solljus	sv	Distributör
tr	Güneş ışığından koruyun	tr	Distribütör
zh	避免阳光直射	zh	分销商

**EC REP**

en	Authorized representative in the European Community
ar	ممثل معتمد في الجمعية الأوروبية
bg	Упълномощен представител в Европейската общност
bs	Ovlašteni predstavnik za Evropsku zajednicu
cs	Autorizovaný zástupce v Evropském společenství
da	Autoriseret repræsentant i Det Europæiske Fællesskab
de	Bevollmächtigter in der Europäischen Gemeinschaft
el	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
es	Representante autorizado en la Unión Europea
fi	Valtuutettu edustaja Euroopan yhteisössä
fr	Mandataire autorisé dans la Communauté européenne
he	הנציג המורשה בקהילה האירופית
hr	Ovlašteni predstavnik u Europskoj zajednici
hu	Meghatalmazott képviselő az EU-ban
it	Mandatario nella Comunità Europea
ko	유럽 공동체의 공인 대리인
nl	Gemachtigd vertegenwoordiger in de Europese Gemeenschap
no	Autorisert representant i den Europeiske Union
pl	Autoryzowany przedstawiciel w Unii Europejskiej
pt	Representante autorizado na Comunidade Europeia
ro	Reprezentant autorizat în Comunitatea Europeană
ru	Уполномоченный представитель в Европейском сообществе
sk	Autorizovaný zástupca v Európskom spoločenstve
sl	Pooblaščen zastopnik v Evropski skupnosti
sr	Ovlašćeni predstavnik u Evropskoj zajednici
sv	Auktoriserad representant i EU
tr	Avrupa Topluluğu yetkili temsilcisi
zh	欧共体授权代表

en	Importer
ar	المستورد
bg	Вносител
bs	Uvoznik
cs	Dovozce
da	Importør
de	Importeur
el	Εισαγωγέας
es	Importador
fi	Maahantuojaja
fr	Importateur
he	היבואן
hr	Uvoznik
hu	Importőr
it	Importatore
ko	수입업자
nl	Importeur
no	Importør
pl	Importer
pt	Importador
ro	Importator
ru	Импортёр
sk	Dovozca
sl	Uvoznik
sr	Uvoznik
sv	Importör
tr	ithalatçı
zh	进口商

English

Indications for Use

Orbit®soft Infusion Sets are intended for the subcutaneous delivery of fluids and medication, such as insulin, from an external infusion pump.

Contraindications

Orbit®soft Infusion Sets are not intended for intravenous (IV) infusion or the infusion of blood or blood products.

Warnings

- ⚠ **Follow the instructions for use** of Orbit®soft Infusion Set and your infusion pump. Use the Orbit®soft Infusion Set with your infusion pump only after consultation with your physician and/or health care professional.
- ⚠ **Avoid infections.** Do not use the Orbit®soft Infusion Set if the individual packaging is damaged, has already been opened, or if the cannula protector is missing. Do not use the Orbit®soft Infusion Sets after the expiration date on the packaging. If these instructions are not followed, there is a risk of infection. If you have any signs or symptoms of an infusion site infection, immediately replace your infusion set and contact your physician for treatment.
- ⚠ **Single use only!** Do not re-use, re-insert and/or re-sterilize the cannula. Re-use and/or re-insertion may cause damage to the cannula and lead to site irritation, infection, and/or inaccurate medication delivery. Re-sterilization may affect the integrity of the Orbit®soft Infusion Set.
- ⚠ **Do not use Orbit®soft Infusion Set for longer than 72 hours.** Longer use may lead to site irritation, infection, and/or inaccurate delivery of medication. Change your infusion site if it becomes itchy, red, inflamed or if you have an unexplained high blood sugar reading.
- ⚠ **Ensure correct delivery of medication.** Before insertion, check that the tubing allows free flow of medication, following the priming instructions of your infusion pump. This will prevent incorrect dosage caused by air in the fluid path or a clogged tubing. Use a new infusion set if priming is not successful. Never prime the infusion set tubing or attempt to free a blocked line while the tubing is connected to the inserted cannula base. You may accidentally deliver too much medication. Do not re-insert the introducer needle into the soft cannula. Reinsertion may cause tearing or puncturing of the soft cannula, which may result in medication leakage.

Replace the Orbit®soft Infusion Set if the tape becomes loose, if the soft cannula becomes kinked during insertion or if the cannula becomes fully or partially dislodged from the skin. Check frequently to make sure the soft cannula remains firmly in place. It must be completely inserted at all times to receive the full amount of medication.

For insulin therapy, 2–3 hours after insertion, a blood glucose test has to be performed to ensure adequate delivery of insulin. Therefore, do not change the Orbit®soft Infusion Set immediately before bedtime.

Develop a conventional insulin injection or medication replacement plan with your healthcare professional if any part of the infusion system fails or system supplies are depleted.

If these instructions are not followed, there is a risk of incorrect delivery of medication. Seek immediate medical assistance if you have any health concern.

⚠ Avoid bending or breaking the cannula/introducer needle.

A broken needle can remain stuck in the body or remain completely under the skin and may require surgical removal.

⚠ For Healthcare Professionals: Avoid recapping or destroying the introducer needle. Recapping or destroying the introducer needle can lead to accidental sharps injuries which can transmit serious infections. Healthcare professionals must follow the specific directions of local regulations for healthcare professionals using infusion sets, including directions for recapping and disposing. Seek immediate medical assistance if you receive a sharps injury.

Caution (please note)

- Orbit®soft Infusion Sets are designed to work with continuous subcutaneous infusion pumps and pump reservoirs with standard luer lock connection. Follow the instructions for use of your infusion pump and/or reservoir.
- Change the infusion site regularly in accordance with your healthcare professional's recommendations. If this is not done, you are at risk for infusion site infection, issues with medication delivery, insulin resistance, and/or lipohypertrophy.
- Consult your healthcare professional to choose the correct cannula length. Your healthcare professional can provide recommendations regarding the specific type of infusion set that is best matched to your needs to ensure accurate delivery of medication.
- Do not put perfumes, deodorants or insect repellent on the Orbit®soft Infusion Set as this may affect the integrity of the infusion set. Protect the product from sunlight and heat. Store Orbit®soft Infusion Sets in a dry location at room temperature and avoid high humidity.

English

Disposal

Dispose of the infusion sets according to your healthcare professional's advice and government regulations for needle/sharps disposal.

Figure 1–2

Choose an infusion site away from the waistline, bones, scar tissue, belly button and recent infusion sites. Wash hands thoroughly and clean the infusion site according to the instructions provided by your healthcare professional. Prepare the infusion site by removing excess hair. We recommend to wipe the site in a circular outward motion, using a 70 % isopropyl alcohol pad, and let it air dry. Make sure the infusion site is completely dry before proceeding. The use of 70 % isopropyl alcohol pads supports the optimal adhesion of the tape. The use of other alcohol pads may result in premature tape failure and interruption of medication delivery.

Figure 3

After connection of the tubing to the infusion pump make sure all air is removed by priming the tubing assembly according to your pump manufacturer's priming instructions.

Figure 4–5

Remove tape backing and be careful not to touch adhesive. Carefully remove the cannula protector from the cannula.

Figure 6–7

Stabilize the infusion site and insert the Orbit®soft cannula with an angle of 90°. For insertion of the cannula, you may also use the Orbit®Insertor from Ypsomed. For maximum adhesion, press the tape to your skin and run your fingers over the adhesive for some seconds.

Figure 8–10

Hold tape against skin with two fingers of one hand. Gently squeeze tabs and remove the introducer needle carefully with the other hand. Dispose the introducer needle into a safe container unit immediately after injection to avoid accidentally pricking.

⚠ You may use the blue needle cover provided to cover the introducer needle after removal (not for healthcare professionals).

**Figure 11 – 13**

Attach the tubing cap straight to the cannula base. One to two clicks should be heard.

△ Rotate the tubing left and right, at least one full turn in each direction, while pulling upward on cap to ensure the tubing cap is fully engaged and the fluid path is opened.

After proper connection of the tubing to the cannula base, the newly inserted cannula has to be filled with medication by priming the cannula head with the priming volumes indicated below. Failure to do so results in missed medication.

Figure A: Disconnection

To disconnect the tubing, gently squeeze tabs on both sides of the tubing cap and remove it. There is no need to cover the infusion base as the septum provides closure.

Before you reconnect the tubing cap, remove the air bubbles from the tubing according to the priming instructions of your infusion pump and swab septum with a 70 % isopropyl alcohol pad. Proper swabbing of the septum reduces germs ingress into the fluid path. To reconnect, follow instructions in figures 11 – 13.

Monitor your blood glucose levels when disconnected and about 2 – 3 hours after reconnecting.

**Priming volume for the tubing assembly
(volumes are approximate):**

45 cm (18 inch): 11 units* (0.11 ml)

60 cm (24 inch): 14 units* (0.14 ml)

75 cm (30 inch): 17 units* (0.17 ml)

105 cm (42 inch): 25 units* (0.25 ml)

**Priming volume for the cannula part (head)
(volumes are approximate):**

6 mm cannula with base: 0.3* units (0.003 ml)

9 mm cannula with base: 0.4* units (0.004 ml)

* U-100 Insulin

Orbit[®]soft products are available as a whole infusion set and as cannula (head) only.

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- SV** Ypsomed AB
- UK** Ypsomed Ltd.
- US** Infirst Healthcare Inc.
Ypsomed AG
Ypsomed Distribution GmbH