
Porelast®Panelast®Fortelast® Plus

Kurzzügige Pflasterbinde
Bande de pansement à allongement court
Short-stretch Adhesive Bandage
Venda adhesiva de tração corta
Ligadura adesiva de baixa elasticidade
Benda adesiva ad allungamento limitato
Korte-rekpleisterzwachtel
Kortfristet plasterforbinding
Lägelastisk klisterbinda
Adhezivní obinadlo s krátkým tahom
Náplastové ovinadlo s krátkým ťahom
Rövid megnyúlású öntapadás pólya
Bandaż przylepny o krótkim naciągu
Kısa çekişli flaster bandaj

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Intended purpose:
For compression, to immobilise limbs, for support and relief

Indications:
Phlebological and lymphological indications: varicosis, venous insufficiency, thrombophlebitis, phlebothrombosis, post-thrombotic syndrome, stage I to III b chronic venous insufficiency according to Widmer, venous and mixed ulcers, oedema (phleboedema, lymph-oedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery. Contusions, sprains, tendon injuries and joint dislocations, healed fractures.

Contraindications:
Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:
In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions: - Incorrect bandaging can result in skin necrosis and pressure damage to nerves. - Slight bluish tinge to the toes is normal, this should disappear as soon as the patient starts to walk. If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression. - The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage. - For use with ABPIs over 0.8. When the ABPI is between 0.5 and 0.8, apply under medical supervision only. - Known allergy and / or hypersensitivity to any of the product components.	
--	--

General instructions: - Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status. - It is recommended to apply the	
--	--

General instructions:
In case of a serious incident, contact the manufacturer and the responsible health authorities.

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