
Raucodrape®

Inzisionsfolie
Champ à inciser
Surgical Drape (Incise Drape)
Lámina de incisión
Película para incisiones
Foglio di incisione
Chirurgisch incisiefolie
Inzisionsfolie
Inzisionsfilm
Incizní fólie
Incizna fólia
Műtési fólia
Folia chirurgiczna

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en Instructions for Use

Product description and performance characteristics

The product is a conformable, self-adhesive surgical incise drape. Natural rubber latex is not an intentional ingredient in the production of the above mentioned product or its packaging. It cannot completely exclude any potential contaminations with latex traces from the production or environment.

Product composition

The product is composed of coated polyurethane film (polyurethane film, adhesive), release paper, adhesive tape.

Intended purpose

The product is intended to cover the operating field during surgical procedures.

Indications

Not applicable.

Contraindications

Do not use the product on patients with known sensitivity to acrylic.

Side effects

In (very) rare cases skin irritations and / or allergies may occur.

Warnings and precautions

During minimally invasive surgical procedures care must be taken that after application of the product introducing parts of the film into the body cavity is avoided by using special instruments.

Application instructions

1. Preparation

- Disinfect the surgical or operating area according to hospital policy or physician-prescribed routine. Allow disinfectant to dry completely. Drape the patient. Affix the product to dry skin only.
- 2. Application**
 - 2.1 Open peel pouch. Remove film sterily.
 - 2.2 Unfold the product and grip the holding edge. It is recommended that the application be performed by two people.
 - 2.3 Peel off cover paper parallel to the holding edge.
 - 2.4 Hold the product a few centimeters above the skin.
- Using a sterile towel, press down firmly on the film, first contacting the skin along the intended incision line to ensure good skin contact. Firmly smooth film into place without tension, working away from the intended incision line to achieve wrinkle free adhesion.
- 2.5 To remove, overstretch the product by gently pulling on the holding edge and remove from the patient's skin.
- 2.6 Do not peel off the product from the skin at a 90o angle.

Storage and transport conditions

Keep dry and protect from light and dust.
The product has a shelf life of 5 years from the date of manufacture.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 1801 and packaging waste a waste code from Chapter 1501 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General Information

In case of serious incident, contact the manufacturer and the responsible health authorities.

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Shanghai ISO Medical Products Co., Ltd.,
Lotus NL B.V.
Lohmann & Rauscher International GmbH & Co. KG