AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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 english

Dear Customer,

thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read and observe these **instructions for use** carefully. If you have any questions, please contact your doctor or medical retailer.

Intended purpose

SofTec Lumbo is a medical device. It is an orthosis¹ for immobilizing the lumbar spine with mobilization function with removable shell.

Stage 1 (Immobilization phase):

Securing the spine after surgery or at the start of non-surgical treatment by using a shell.

Stage 2 (Mobilization phase):

To allow for mobilization over the course of recovery, the shell is removed and replaced by CFRP rods (carbon fiber reinforced plastic rods) in the knitted fabric.

Stage 3 (Activity phase):

With increasing activity, the CFRP rods are replaced by corset stays, thereby reducing the external supportive effect.

 **SofTec Lumbo must be individually adjusted for the patient by a trained specialist². This is the only way the full effectiveness¹ and ideal wearing comfort of the orthosis are guaranteed. If the product is not used properly any claim for warranty will expire.**

Indications

- Spondylolysis / spondylolisthesis, grades II – III

- Extremely severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions/herniation
- Most severe (pseudo) radicular lumbar syndrome
- Degenerative facet syndrome
- Condition after herniation, e.g. post-operative care
- Condition after disectomy
- Condition after spondylosis / kyphoplasty (one or more levels)
- Lumbar spinal canal stenosis
- Vertebral fractures (lumbar spine) with vertebral body damage
- Tumors (metastases)
- Severe degeneration / muscular insufficiency

Risks of using this product

 **Important information**

Please closely observe the specifications in these instructions for use and the instructions given by the specialist².

- SofTec Lumbo is a prescribed product that should be worn under a physician's guidance. SofTec Lumbo should only be worn in accordance with the instructions contained in these usage directions and for the listed areas of application.
- The plastic shell or stays of the SofTec Lumbo may only be modified by a properly trained specialist. Failure to comply with this requirement may adversely affect product performance, thereby canceling any product liability.
- If the knitted sections of SofTec Lumbo are subsequently modified, the proper functioning of the support can no longer be guaranteed.
- Should you notice any unusual changes (e.g. increased symptoms), please contact your doctor without delay.
- Please consult your doctor before using SofTec Lumbo in combination with other products.
- Do not allow the product to come into contact with ointments, lotions or substances containing grease or acids.
- SofTec Lumbo can readily be disposed of with domestic refuse.
- Side effects, with the assumption of correct fit involving the body as a whole have not been reported to date. Any supports or orthoses applied externally to the body can, if tightened excessively, signs of localized pressure, rarely, constrict the underlying blood vessels or nerves. In such cases, please loosen the straps or the fastening slightly and have the correct fit and size of the orthosis checked if necessary.
- Caution: When further loosening the straps or tensioning and fastening systems, or when removing the product or adjusting configured limitations of movement, there is not sufficient support / sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.

 Contains natural rubber latex (elastodiene), which can cause skin irritation in persons with an allergy. Take off the orthose if you notice any allergic reaction.

Contraindications

- Hypersensitive reactions harmful to health have not been reported to date. In the following conditions such supports or orthoses should only be fitted and worn after consultation with your doctor:
- Skin disorders /injuries to relevant parts of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness and excessive heat build-up.
 - Impaired sensation.
 - Severe restriction in the performance of heart and lungs (risk of elevated blood pressure if the aid is worn during severe physical exertion).

Application instructions

Fitting and removing SofTec Lumbo

Despite its special design, SofTec Lumbo is always fitted and removed in the same way at all stages.

Putting on the orthosis:

- The shell individually adapted by a specialist² must sit loosely on the iliac crests in the pelvic area (Picture 2). When putting on the orthosis without the shell, it must be ensured that the support is positioned correctly. If necessary, practice correct positioning with the treating specialist.
- Place the support around the trunk so that the reinforcements visible on the outer knitted section are positioned parallel to the spine (Fig. 3).
- Place your fingers in the specially provided loops on the two halves of the light grey fastener and pull evenly forward (Fig. 5).

- Press the left half of the fastener against your stomach and place the right half above it to secure the Velcro fastening. If the optional pressure pad for the abdominal wall is used, ensure that this is positioned correctly (this pressure pad should be centrally positioned beneath the fastening).
- When securing the Velcro fastening, remove your fingers from the left loop first, followed by the right loop (Fig. 6).
- Take hold of the two LPT (**Low Power Tension**) straps in both hands and pull evenly forward (Fig. 7). Secure both LPT straps onto the large fastening (Fig. 8).

- With the "long" version, always fasten the lower straps first and then secure the upper straps on top.
- During the first fitting of the orthosis, the optimal, individual fit must be checked, and corrected if necessary, by a trained specialist. Correct fitting should be practiced.

Removal:

- Undo the Velcro fastenings of the LPT straps on the outer knitted section and pull back the straps to loosen the support and fasten the ends of the straps to the retracted sections.
- Open the large light grey fastener.
- Remove the support and reclose the light grey fastener.

Product care

 **Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!**

- SofTec Lumbo can be washed separately by hand at 30 °C (86 °F) using a mild detergent. Ensure that the plastic shell is removed beforehand. Proceed as follows:
- Undo the Velcro fastenings of the two thin elastic straps. Pull the straps out of the holes and remove.
- Undo the two Velcro fastenings between the outer and inner knitted sections located at the center of the back.
- Remove the plastic shell from the support.
- The shell can easily be washed with a sponge, using a pH-neutral soap or shampoo. Do not dry using direct heat (dryer, heater). Either rub dry with a cloth or air dry.
- Removal of the CFR or corset rods is not necessary.
- If the rods are to be removed from the outer knitted section, turn back the fasteners at the top (red section) and push out the rods (Fig. 9).
- To avoid damage and to preserve their long-term adhesive strength, we recommend closing the Velcro fastenings before washing (including the small Velcro dots between the inner and outer knitted sections).
- After washing, perform the above sequence in reverse. When re-inserting the CFR rods, check the markings for correct orientation (top / bottom and left / right). Also check the correct positioning of the knitted sections containing the shell. Fasten identically coloured Velcro dots to each other (the sewn-in labels are always at the top, the shell is marked towards the top).

This product has been examined by our own quality control system. However, should you have any complaint, please contact your retailer or our customer service department.

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SofTec Lumbo have not been observed. Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a medical specialist. Before first using our medical product, please seek the advice of a physician or trained specialist, as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this specialist as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the specialist, please contact your physician or retailer, or contact us directly.

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Polyethylene (PE), Polyamide (PA), Polyurethane (PUR), High Grade Steel, Elastodiene (ED), Ethylene propylene diene monomer rubber (EPDM), Cotton (CO), Polyoxymethylene (POM), Spandex (EL), Polyester (PES), Glass fiber (GF), Thermoplastic polyurethane (TPU), Carbon (CF)

 – Medical Device

 – DataMatrix Unique Device Identifier UDI

Version: 2023-02

¹Orthosis = orthopaedic appliance used to stabilize, relieve, immobilize, control or correct the limbs or the torso

²A specialist is any person who is authorized according to the state regulations for fitting and instruction in the use of supports that are relevant to you.

Specialists²:

Application instructions

Fitting and removing SofTec Lumbo

Despite its special design, SofTec Lumbo is always fitted and removed in the same way at all stages.

Fitting:

- Wrap the orthosis around your body so that the inserts attached to the outer knitted fabric for the CFRP and corset stays are positioned parallel to the spine (Fig. 1), sewn-in label at the top.
- The anatomical individually adapted shell must be positioned by the specialist evenly above the iliac crests (Fig. 2).
- When the lower knitted section is being closed (if possible by the patient) ensure that the support is correctly centered (Fig.3).
- The individually adjusted shell with the precontoured pelvic indentations must must sit loosely on the iliac crests. When putting on the orthosis without the shell, it must be ensured that the support is positioned correctly (label at the top).
- For individual adaptation of the shell, the wedges can be cut into or cut out and the shell adapted accordingly (Fig. 4).
- Place the fingers in the specially provided loops on the two halves of the light grey fastener and pull evenly forward (Fig. 5).
- Press the left half of the fastener against the patient's stomach and pull the right half over it to secure the Velcro fastening (Fig. 6). Ensure that the pressure pad for the abdominal wall, if inserted, is positioned correctly (this pressure pad should be fastened in the center of the left, underlying Velcro fastening).
- When securing the Velcro fastening, remove the fingers from the left loop first, followed by the right loop.
- Take hold of the two LPT (**Low Power Tension**) straps in both hands and pull evenly forwards until the correct level of tension is reached (Fig. 7). Secure both LPT straps onto the large fastening (Fig. 8). With the "long" version, always fasten the distal (lower) straps first, followed by the proximal (upper) straps. The plastic buckles on the LPT straps can be shaped by thermoplastic formation (i.e. with the application of heat) if necessary.
- Check, and if necessary correct, the optimal, individual fit.
- The CFRP stays in the outer knitted fabric are thermoformable (using heat). To remove the rods (also when changing the corset rods), open the inserts at the top. When re-inserting the CFR rods, check the markings for correct orientation (top / bottom and left / right) (Fig. 9).
- The patient should practice the correct way to put on the orthosis with help from the specialist.

Removal:

- Undo the Velcro fastenings of the LPT straps on the outer knitted section and pull back the straps to loosen the support and fasten the ends of the straps to the retracted sections.
- Undo the two halves of the light grey fastener.
- Remove the support and reclose the light grey fastener.

Adapting SofTec Lumbo

- The shell must be adapted to the patient's individual anatomy. Pressure points must be avoided at all costs. The plastic shell can be heated to approx. 80 °C (176 °F) using a hot-air gun.
- The CFR rods can also be heated and reshaped.

 **Note: When SofTec Lumbo is used after stabilizing operations (spondylolyses), the CFR rods must be shaped to match the normal spinal curvature (angle of anatomical lordosis) in order to avoid the transfer of excessive straightening forces (lordosis flattening).**

- When using SofTec Lumbo for bridging the lordosis, insert the straight CFR rods as supplied with the convex side facing outwards (instead of shaping them to match the body's contours as for stabilization of the lordosis).
- If bending hooks are employed, only round bending hooks should be used.
- Instructions for removing the shell and rods are provided under "Product care".
- The ideal temperature for adapting the stays is between 170 °C and 190 °C. Heating above these temperatures leads to destruction of the material.
- When heating carbon / glass-fibres keep the material a minimum distance of 1,97 inch from the heat source to avoid burning the surface.
- The carbon / glass-fibre should be heated evenly across the surface to avoid damaging the material.

 **In that case Bauerfeind will not be liable for the product.**

- The material should only be bent gently and in large arches.
- The cooling process can be accelerated using cold water which will not change the material's characteristics.
- The material can be heated and molded repeatedly.
- The heating process does not visually alter the surface of the material. In order to find out the appropriate time for bending, test it during the heating process.
- Check the fit of the support on the standing and sitting patient after each modification.
- Local pressure points exerted on the skin or around the iliac crests must be avoided at all costs.

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