



Zahnseidenhalter_169500001

Zahnseide_494099070



So einfach geht's:

Geöffneten Halter in die Aussparungen der Spezialbox legen (Seide wird durch Schließen des Verschlusshebels gespannt), an der Schneidvorrichtung abtrennen, fertig. Nach jedem Zahnzwischenraum frische Zahnseide verwenden.

Here's how it works:

Simply insert the open holder into the slots on the special box (the floss is tightened by closing the lever) and detach using the cutter. That's it! Use a fresh section of floss between each pair of teeth.

(1) Das Set enthält den **ergonomisch geformten Cleudent-Halter** und **eine Box mit Qualitäts-Zahnseide**. Das System funktioniert natürlich auch mit allen handelsüblichen Zahnseiden.

(2) Einzigartig **einfaches Spannsystem** für perfekte Reinigungswirkung

Alle Zahnzwischenräume – auch die der Backenzähne – werden mühelos erreicht.

Hoher Hygienestandard – keine Berührung der Zahnseide

Größte Effizienz bei geringem Verbrauch

Für jedes Alter geeignet (Empfehlung ab 10 Jahren)

(1) The set includes the **ergonomically-shaped Cleudent holder** and a **box with high-quality dental floss**. Naturally the system also works with all commercially obtainable dental floss.

(2) Unique and **simple tensioning system** for perfect cleaning efficiency

You can effortlessly reach all interdental spaces – even between the molars.

High hygienic standard – no contact with floss

Maximum efficiency and low consumption

Suitable for all ages (recommended from 10 years)

Farbe

Colour

