

Melgisorb® Ag

Highly absorbent
antimicrobial alginate
dressing

Silver-containing alginate

- Provides sustained antimicrobial effect
- In vitro testing show antimicrobial effect from 4 hours up to 14 days



Alginate with CMC fibers

- Suitable for high exuding wounds
- Designed for shallow and deep wounds
- Allows removal in one piece

- Excellent absorption for heavily exuding deep wounds
- Provides sustained antimicrobial effect
- High fibre strength allows removal in one piece

How to use Melgisorb® Ag



1. Debride and irrigate the wound site in accordance with standard protocols. Select a size of Melgisorb® Ag that is slightly larger than the wound. Cut or fold the dressing to fit the wound.



2. Loosely fill deep wounds, ensuring the dressing does not overlap the wound margins. Apply to wound bed directly.



3. Cover and secure Melgisorb® Ag with a secondary dressing.

Areas of use

Melgisorb® Ag is indicated for the management of moderate to heavily exuding partial to full thickness wounds such as:

- Post-operative surgical wounds
- Trauma wounds (dermal lesions, trauma injuries or incisions)
- Leg ulcers
- Pressure ulcers
- Diabetic ulcers
- Graft and donor sites
- Partial thickness burns

This product can also be used under compression bandages. As the product contains alginate, it may assist in supporting the control of minor bleeding in superficial wounds. Melgisorb® Ag may be used in the management of infected wounds or wounds in which there is an increased risk of infection (under the medical supervision of a healthcare professional). Melgisorb® Ag is indicated for external use only.

Note

Dressing change frequency will depend on the wound condition and the level of exudate. Initially, it may be necessary to change the dressing every 24 hours. If the wound appears dry, saturate the dressing with sterile saline solution prior to removal. Gently remove the dressing from the wound bed. Irrigate the wound site in accordance with standard protocols prior to application of a new dressing.

Contra indications

Melgisorb® Ag is not indicated for use on the following:

- Individuals with a known sensitivity to alginates or silver
- To control heavy bleeding
- For direct application on dry or lightly exuding wounds.

For detailed information, please see instructions for use.

Assortment (sterile packed)

Art. No.	Size (cm)	Pack count box/trp
256055	5 x 5	10/100
256105	10 x 10	10/100
256155	15 x 15	10/100
256205	20 x 30	5/50
256605	3 x 44	10/100



Find out more at molnlycke.com

Mölnlycke Health Care AB, P.O. Box 6, Entreprenörstråket 21, SE-431 21 Mölndal Sweden. Phone + 46 31 722 30 00. The Mölnlycke and Melgisorb trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. ©2025 Mölnlycke Health Care AB. All rights reserved. HOIM005898.

